



# INNOVATION IRRIGATION™ SYSTEMS


**ATC TECHNOLOGIES, Inc.**, the market leader in Suction Irrigation Products, brings the best selection of Cannulas, Electrosurgical Probes, and the most comprehensive Accessories in Fluid Management.

Select from Pediatrics in 21cm, to General Laparoscopy in 28cm & 33cm, to Obesity Procedures using 45cm lengths

## Back to Basics...


### 3mm Cannula Probe Tip

For optimum M.I.S. procedures. With or without irrigation holes.




### 5mm Cannula Probe Tip

The standard in suction irrigation. With or without irrigation holes.



### 5mm Atraumatic Probe Tip

Allows for delicate aspiration. With or without irrigation holes.



### 10mm Cannula Probe Tip

Largest lumen for optimum aspiration. With or without irrigation holes.




### 5mm Composite Probe Tip

Non-conductive for use with Monopolar electrode inserts. With or without irrigation holes.



### 10mm Composite Probe Tip


Non-conductive for use with Monopolar electrode inserts. With or without irrigation holes.



## Specialties...


### 5mm Cannula Micro Tip

A reduced tip designed to provide higher pressure for dissection.




### 5mm Cannula Pool Tip

Allows for Atraumatic aspiration of fluids without suctioning loose or delicate tissue.




### 5mm Hydro-Scout™ Probe Tip

Duckbill curved tip maneuvers well around tissue, organs, and vessels for dissecting and debridement.



### 10mm Cannula Pool Tip

Large volume Atraumatic aspiration of fluids without suctioning loose and delicate tissue.




### 10mm Stone Scoop Probe Tip

Unique scoop design to retrieve stones and specimens. Drainage holes and sheath allow capsulation and cleaning of specimens.




### 10mm Instrument Insert Probe Tip

Non-conductive for use with Electrosurgical instruments, electrodes, and accessories. With or without irrigation holes.



## Innovation...

Threaded base for Tight seal and no wobble.



## Innovation...

Universal Seals for pass thru instruments, electrodes, and accessories.

