


# Four-in-one Electrocautery Probes




## 33cm Length x 5mm Diameter Disposable Electrosurgical Probes

	Reorder #	Product Description	Quantity
	4180	Spatula Tip Style Electrode with Sliding Sheath	10/box   Sterile
	4181	J-Hook Tip Style Electrode with Sliding Sheath	10/box   Sterile
	4182	L-Hook Tip Style Electrode with Sliding Sheath	10/box   Sterile
	4183	Needle Tip Style Electrode with Sliding Sheath	10/box   Sterile
	4184	Ball Tip Style Electrode with Sliding Sheath	10/box   Sterile
	4185	Blunt Tip Style Electrode with Sliding Sheath	10/box   Sterile

## 45cm Length x 5mm Diameter Disposable Electrosurgical Probes

	Reorder #	Product Description	Quantity
	4186	Spatula Tip Style Electrode with Sliding Sheath	5/box   Sterile
	4187	J-Hook Tip Style Electrode with Sliding Sheath	5/box   Sterile
	4188	L-Hook Tip Style Electrode with Sliding Sheath	5/box   Sterile
	4189	Needle Tip Style Electrode with Sliding Sheath	5/box   Sterile
	4190	Ball Tip Style Electrode with Sliding Sheath	5/box   Sterile
	4191	Blunt Tip Style Electrode with Sliding Sheath	5/box   Sterile

## Reusable Accessories

	Reorder #	Product Description	Quantity
	4103	10 ft. Monopolar High Frequency Cord, Female End to Bovie End	1/box   Non-Sterile
	4104	10 ft. Monopolar High Frequency Cord, 90° Female End to Bovie End	1/box   Non-Sterile
	4199	Closure Plug	2/box   Non-Sterile



30-B Upton Drive | Wilmington, MA 01887  
 781.939.0725 | Fax: 781.939.0726  
[www.atctechnologiesinc.com](http://www.atctechnologiesinc.com)