


CROSSBOW™

Suture ease

Premium Suture Passer

- ✓ **Easy Loading**
 - Large aperture
- ✓ **Secure Holding**
 - Braided or monofilament suture
- ✓ **Positive Release**
 - Exclusive suture ejection feature



Versatile Guide

Robust Mechanism

- Intuitive - one moving part
- Robust "feel"

User Friendly Design

- Sleek guide tip for smooth entry
- Low insertion force

Visual Indicator


- Reliable guide placement

Snare Tip Landmark

- Enables predictable suture placement
- Consistent fascial purchase

Extra Large Snare Loops

- Fast & easy suture retrieval
- Reliable target zone
- Very strong nylon



Product Code

Description

SE-CB

CrossBow Fascial Closure System,
5 Pack Dispenser Box