



CROSSBOW™ Fascial Closure System


Simple Suture Retrieval



Place first tail of suture through 1st snare




Place other tail of same suture through 2nd snare



Retract snares to retrieve suture tails

Reproducible Suture Placement



Visual indicators allow for accurate positioning of the guide relative to the fascia, ensuring precise suture placement

Adaptor (optional)

Adaptor (optional)
Designed for 15mm and larger trocar defects, offering versatility of multiple suture patterns

- Ideal for specimen removal sites
- 18.5 - 21mm Diameter
 - 105 mm working length
 - Simple "snap on" - "twist off"

EzStitch™ Port Site Closure System

Exceptional Safety

Needle protected within distal chamber of guide

Simple Suture Retrieval


Suture is automatically captured in distal chamber

Secure Fascial Closure

Reliable and efficient suture placement

Notch

For repeatable fascial edge placement



Sharps Protection

Easy load
Suture Passer

Adjustable
Sealing Cone
accommodates larger defects

Extended length
Guide Shaft
ideal for bariatric applications

Distal Chamber
with auto suture capture





Product Code Description

- SE-CB** CrossBow™ Fascial Closure System
5 Pack Dispenser Box
- SE-CB-XLA** CrossBow™ Fascial Closure System w/ ADAPTORTM
5 Pack Dispenser Box

SecurusEP

Product Code Description

- SE-EP** SecurusEP Suture Passer
10 Pack Dispenser Box




EzStitch™

Port Site Closure System

Product Code Description


- SE-EZS** EzStitch Port Site Closure System
5 Pack Dispenser Box
- SE-EZDG** EzStitch Disposable Guide
10 Pack Dispenser Box
- SE-EZRP** EzStitch Reusable Suture Passer
Individual



Suture Ease

... when you need closure.


Simplify your laparoscopic port site closure technique



INNOVATION OF THE YEAR



Fast
Easy
Versatile



©2017 Suture Ease, Inc. All rights reserved. CrossBow Fascial Closure System is a trademark of Suture Ease, Inc.

Suture Ease, Inc.
1735 N. First St.,
Ste 300
San Jose, CA 95112

www.sutureease.com
info@sutureease.com

Customer Service
408-459-7595